# The LIFETIMES Study – Interim Results from a Health and Economic Evaluation of a Rapid Phenotypic AST System


Tiziana Di Martino<sup>1</sup>, Gian Maria Rossolini<sup>2</sup>, Fabio Gori<sup>4</sup>, Fabio Sanguinetti<sup>3</sup>, Antonella Repetto<sup>4</sup>, Alberto Antonelli<sup>2</sup>, Giulia De Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Gori<sup>4</sup>, Paola Bondi<sup>4</sup>, Fabio Morecchiato<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Gori<sup>4</sup>, Paola Bondi<sup>4</sup>, Fabio Morecchiato<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Gori<sup>4</sup>, Paola Bondi<sup>4</sup>, Fabio Morecchiato<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Gori<sup>4</sup>, Paola Bondi<sup>4</sup>, Fabio Morecchiato<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Morecchiato<sup>3</sup>, Chiara Chilleri<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Gori<sup>4</sup>, Paola Bondi<sup>4</sup>, Fabio Morecchiato<sup>2</sup>, Chiara Chilleri<sup>2</sup>, Journa Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>2</sup>, Bruno Viaggi<sup>2</sup>, Fabio Songuinetti<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Giani<sup>3</sup>, Angelo Giacomucci<sup>4</sup>, Tommaso Gi Brunella Posteraro<sup>3</sup>, Giorgia Sanna<sup>3</sup>, Fabiola Cammarota<sup>3</sup>, Giulia Menchineli<sup>3</sup>, Filippo Allegrucci<sup>4</sup>, Riccardo Accattoli<sup>4</sup>, Venere Cortazzo<sup>5</sup>, Marilena Agosta<sup>5</sup>, Francesca Ippolita Calò Carducci<sup>5</sup>, Anna Bakidou<sup>1</sup>, Leonie Vogt<sup>1</sup>, Eric Eriksson<sup>1</sup>, Martin Svedberg<sup>1</sup>, Oliver Mainwaring<sup>1</sup>, Marco Barbieri<sup>6</sup>, Michael Drummond<sup>6</sup> <sup>1</sup>Q-linea, Uppsala, Sweden.<sup>2</sup>Azienda Ospedale Pediatrico Bambino Gesù IRCCS, Rome, Italy. <sup>6</sup>Centre for Health Economics, University of York, York, United Kingdom <sup>1</sup>Q-linea, Uppsala, Sweden.<sup>2</sup>Azienda Ospedale Santa Maria della Misericordia, Perugia, Italy. <sup>4</sup>Azienda Ospedale Pediatrico Bambino Gesù IRCCS, Rome, Italy. <sup>4</sup>Azienda Ospedale Pediatrico Bambino Gesù IRCCS, Rome, Italy. <sup>4</sup>Azienda Ospedale Santa Maria della Misericordia, Perugia, Italy. <sup>4</sup>Azienda Ospedale Pediatrico Bambino Gesù IRCCS, Rome, Italy. <sup>4</sup>Azienda Ospedale Santa Maria della Misericordia, Perugia, Italy. <sup>4</sup>Azienda Ospedale Santa Misericordia, Perugia, Italy. <sup>4</sup>Azienda Ospedale S

### Introduction

The fully automated ASTar<sup>®</sup> System tests directly from positive blood cultures and delivers AST results in approximately six hours<sup>1</sup> (Figure 1).

The LIFETIMES Health Economics and Outcomes Research (HEOR<sup>2</sup>) study is investigating the clinical impact of ASTar when managing ICU patients with Gram-negative Bloodstream Infections (BSIs) at four Italian sites with high incidence of antimicrobial resistance. Here, we present interim data from the LIFETIMES study.





- Clinicians act upon ASTar results by changing the antibiotic course or confirming appropriate empiric
- therapy. ASTar has been compared with standard of care and routine automated AST methods.
- So far, 48 patients have been enrolled 48 were used for performance evaluation and 37 for clinical impact evaluation. Patient characteristics are shown in Table 1.

Table 1. Patient characteristics for all individuals in the clinical impact evaluation.

Patient characteristics		
Total patients (n)		
n = 37		
Age		
Median (IQR)	62.5 (46.8 to 73)	
Minimum–maximum	0—85	
Gender		
Male	29 (78%)	
Female	8 (22%)	
Score (median)		
SAPS (IQR)	40 (31.5 to 54)	
SOFA (IQR)	8 (5.3 to 10)	
CCI (IQR)	3 (2 to 5.8)	
Time in hospital (days)		
Median (IQR)	20 (12 to 45)	
Minimum–maximum	1—78	
Time in ICU (days)		
Median (IQR)	17.5 (8.5 to 34)	
Minimum–maximum	1–49	


### www.qlinea.com

Q-linea and ASTar are registered trademarks of Q-linea AB. Patents: www.qlinea.com/patents

© Q-linea AB 2024. Q-linea AB, Dag Hammarskjölds väg 52 A, SE 752 37 Uppsala, Sweden.


Products are CE-IVD marked and FDA 510(k) cleared\*. Availability of product in each country depends on local regulatory marketing authorization status. \*Data presented above is generated with CE-IVD marked version of ASTar BC G– Kit software. Panel and performance differs from the FDA cleared version. D85636 2024-05-31





# ASTar delivers actionable results >30 hours faster than BMD


Patient samples underwent AST using ASTar and traditional BMD. Median time from +BC to ASTar results was 11.3 h (6.3–23.2 h). Median time from +BC to BMD results was 41.3 h (25.1–82.0 h). The overall median time difference between the methods was 30.0 h (Fig 4).



### ASTar drove antibiotic change in 43% of patient cases

For the 37 patients with complete case reports, we investigated how timely AST results via ASTar influenced treatment adjustment. In 43% (16/37) of patients, ASTar prompted antibiotic changes, with the type of antibiotic change dependent on the individual case.

As demonstrated, ASTar can accurately guide multiple types of antibiotic adjustment. In 57% (21/37) of patients, ASTar confirmed prior therapy (Fig 5).



### ASTar performs at 94% total EA and above 94% total CA

For all tested antibiotics, MIC data was interpreted and assessed using EUCAST breakpoints v13.0<sup>3</sup>.

ASTar maintained a high agreement with reference BMD methods, as assessed by an overall EA of 94.0%, CA of 94.3%, VMD rate of 3.1%, and MD rate of 1.9% (Table 2 and 3).

**Table 2.** Overall performance data. Essential Agreement
 (EA), Categorical Agreement (CA), Very Major Discrepancy (VMD), and Major Discrepancy (MD).

ASTar vs. traditional BMD				
EA	CA	VMD	MD	
#/tot (%)	#/tot (%)	#/tot (%)	#/tot (%)	
584/621	578/613	3/98	9/485	
(94.0%)	(94.3%)	(3.1%)	(1.9%)	

Note: Performance data above is generated with the CE-IVD marked version of ASTar BC G– Kit software. Panel and performance differs in the FDA cleared version of the ASTar BC G-Kit software available in US.

### **Table 3.** Performance data for all antibiotics used in the study.

Antimicrobial agent	ASTar vs. trae	ditional BMD
	EA #/tot (%)	CA #/tot (%)
Amikacin	44/47 (93.6%)	45/47 (95.7%)
Amoxicillin-Clavulanic acid	20/23 (87%)	20/23 (87%)
Ampicillin	3/3 (100%)	3/3 (100%)
Aztreonam	12/12 (100%)	12/12 (100%)
Cefepime	37/39 (94.9%)	36/39 (92.3%)
Cefotaxime	9/9 (100%)	9/9 (100%)
Ceftazidime	32/41 (78%)	31/41 (75.6%)
Ceftazidime-Avibactam	41/44 (93.2%)	43/44 (97.7%)
Ceftolozane-Tazobactam	42/43 (97.7%)	43/43 (100%)
Ceftriaxone	24/26 (92.3%)	26/26 (100%)
Cefuroxime	5/5 (100%)	5/5 (100%)
Ciprofloxacin	44/47 (93.6%)	44/47 (93.6%)
Colistin	17/19 (89.5%)	19/19 (100%)
Ertapenem	31/31 (100%)	31/31 (100%)
Gentamicin	34/36 (94.4%)	35/36 (97.2%)
Levofloxacin	34/34 (100%)	23/26 (88.5%)
Meropenem	45/48 (93.8%)	45/48 (93.8%)
Piperacillin-Tazobactam	42/44 (95.5%)	40/44 (90.9%)
Tigecycline	7/7 (100%)	7/7 (100%)
Tobramycin	22/24 (91.7%)	22/24 (91.7%)
Trimethoprim-Sulfamethoxazole	39/39 (100%)	39/39 (100%)

## Conclusion

- ASTar had an impact on 16 out of 37 patients (43%)
- ASTar shortens the median time from +BC to AST result by 30.0 hours when compared to traditional BMD and 16.3 hours when compared to other automated systems
- ASTar performance aligns with traditional BMD overall EA and CA >90%

### References

- Göransson, J et al. "Performance of a System for Rapid Phenotypic Antimicrobial Susceptibility Testing of Gram-Negative Bacteria Directly from Positive Blood Culture Bottles." Journal of clinical microbiology vol. 61,3 (2023): e0152522. doi:10.1128/jcm.01525-22
- 2. ISPOR. https://www.ispor.org/heor-resources/about-heor; 2024
- . EUCAST. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters, version 13.0, 2023.