

Reportable ranges for antimicrobials and bacteria covered by ASTar BC G– Kit (US). Dots indicate where MIC results and interpretations (when applicable) will be presented.

| Antimicrobial class | Antimicrobial agent | <i>A. baumannii</i> complex | <i>C. freundii</i> complex | <i>C. koseri</i> | <i>E. cloacae</i> complex | <i>E. coli</i> | <i>K. aerogenes</i> | <i>K. oxytoca</i> | <i>K. pneumoniae</i> group | <i>M. morganii</i> | <i>P. mirabilis</i> | <i>P. vulgaris</i> | <i>S. marcescens</i> | <i>P. aeruginosa</i> |
|-----------------------------|--|-----------------------------|----------------------------|------------------|---------------------------|----------------|---------------------|-------------------|----------------------------|--------------------|---------------------|--------------------|----------------------|----------------------|
| Penicillin | Ampicillin | | | | | • | | | | | • | | | |
| β-lactam combination agents | Ampicillin-sulbactam ¹ | • | | • | | • | | • | • | • | • | • | | |
| β-lactam combination agents | Ceftolozane-tazobactam | | • | • | • | • | | • | • | | • | • | • | • |
| β-lactam combination agents | Ceftazidime-avibactam ² | | • | • | • | • | • | • | • | • | • | • | • | • |
| β-lactam combination agents | Meropenem-vaborbactam ³ | | • | • | • | • | • | • | • | • | • | • | • | • |
| β-lactam combination agents | Piperacillin-tazobactam ⁴ | • | • | • | | • | • | | • | • | • | • | • | |
| Cephalosporin | Cefazolin | | | • | | • | | • | • | | • | | | |
| Cephalosporin | Cefepime | | • | • | • | • | • | • | • | | • | • | • | • |
| Cephalosporin | Cefotaxime | | • | • | • | • | | • | • | | • | • | • | |
| Cephalosporin | Ceftriaxone | | • | • | • | • | • | • | • | | • | • | • | |
| Cephalosporin | Cefoxitin | | | • | | • | | • | • | | • | | | |
| Cephalosporin | Cefuroxime | | | • | | • | | • | • | | • | | | |
| Cephalosporin | Ceftazidime | • | • | • | • | • | • | • | • | • | • | • | • | |
| Monobactam | Aztreonam | | • | • | • | • | • | • | • | • | • | • | • | • |
| Carbapenem | Ertapenem | | • | • | • | • | | • | • | | • | • | • | |
| Carbapenem | Meropenem | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Aminoglycoside | Gentamicin | | • | • | | | | • | • | • | • | • | • | |
| Aminoglycoside | Tobramycin | | • | • | • | • | • | • | • | | • | • | | • |
| Aminoglycoside | Amikacin | • | • | • | • | • | • | • | • | • | • | • | • | |
| Tetracycline | Tigecycline | | • | • | • | • | • | • | • | | | | • | |
| Fluoroquinolone | Ciprofloxacin | | • | • | • | • | • | • | • | • | • | • | • | • |
| Fluoroquinolone | Levofloxacin | | • | • | • | • | • | • | • | • | • | • | • | • |
| Miscellaneous | Trimethoprim-sulfamethoxazole ⁵ | | • | • | • | • | • | • | • | • | | • | | |

¹ Ampicillin-sulbactam in the ratio 2:1

² For susceptibility testing purposes, the concentration of avibactam is fixed at 4 µg/mL

³ For susceptibility testing purposes, the concentration of vaborbactam is fixed at 8 µg/mL

⁴ For susceptibility testing purposes, the concentration of tazobactam is fixed at 4 µg/mL

⁵ Trimethoprim:sulfamethoxazole in the ratio 1:19

Antimicrobial reportable ranges for AST and FDA-recognized/approved clinical breakpoints, i.e. susceptible (S), intermediate (I), susceptible-dose dependent (SDD) and resistant (R).

| Antimicrobial agent | ASTar BC G- Reportable range (µg/mL) | Organism Group | FDA-Recognized / Approved Breakpoints (µg/mL) | | |
|-------------------------------|---|-----------------------------|---|------------------|------|
| | | | S | I/SDD | R |
| Amikacin | ≤2 to ≥256 | <i>A. baumannii</i> complex | ≤8 | 16 | ≥32 |
| | ≤0.5 to ≥256 ¹ | Enterobacterales | ≤4 | 8 | ≥16 |
| Ampicillin | ≤1 to ≥128 | Enterobacterales | ≤8 | 16 | ≥32 |
| Ampicillin-sulbactam | ≤1 to ≥128 ² | <i>A. baumannii</i> complex | ≤8 | 16 | 32 |
| | | Enterobacterales | | | |
| Aztreonam | ≤0.25 to ≥128 ³ | Enterobacterales | ≤4 | 8 | ≥16 |
| | ≤0.5 to ≥128 | <i>P. aeruginosa</i> | ≤8 | 16 | ≥32 |
| Cefazolin | ≤0.25 to ≥32 | Enterobacterales | ≤2 | 4 | ≥8 |
| Cefepime | ≤0.25 to ≥128 | Enterobacterales | ≤2 | 4-8 ⁴ | ≥16 |
| | | <i>P. aeruginosa</i> | ≤8 | 16 | ≥32 |
| Cefotaxime | ≤0.25 to ≥16 | Enterobacterales | ≤1 | 2 | ≥4 |
| Cefoxitin | ≤1 to ≥128 | Enterobacterales | ≤4 | 8 | ≥16 |
| Ceftazidime | ≤0.5 to ≥128 | <i>A. baumannii</i> complex | ≤8 | 16 | ≥32 |
| | ≤0.25 to ≥128 ⁵ | Enterobacterales | ≤4 | 8 | ≥16 |
| Ceftazidime-avibactam | ≤0.125 to ≥64 ⁶ | Enterobacterales | ≤8 | - | ≥16 |
| | ≤0.125 to ≥64 | <i>P. aeruginosa</i> | | | |
| Ceftolozane-tazobactam | ≤0.25 to ≥64 | Enterobacterales | ≤2 | 4 | ≥8 |
| | | <i>P. aeruginosa</i> | ≤4 | 8 | ≥16 |
| Ceftriaxone | ≤0.25 to ≥16 | Enterobacterales | ≤1 | 2 | ≥4 |
| Cefuroxime | ≤1 to ≥128 | Enterobacterales | ≤8 | - | ≥16 |
| Ciprofloxacin | ≤0.125 to ≥16 | Enterobacterales | ≤0.25 | 0.5 | ≥1 |
| | | <i>P. aeruginosa</i> | ≤0.5 | 1 | ≥2 |
| Ertapenem | ≤0.06 to ≥16 | Enterobacterales | ≤0.5 | 1 | ≥2 |
| Gentamicin | ≤0.25 to ≥64 | Enterobacterales | ≤2 | 4 | ≥8 |
| | | | | | |
| Levofloxacin | ≤0.125 to ≥32 | Enterobacterales | ≤0.5 | 1 | ≥2 |
| | | <i>P. aeruginosa</i> | ≤1 | 2 | ≥4 |
| Meropenem | ≤0.06 to ≥128 | <i>A. baumannii</i> complex | ≤2 | 4 | ≥8 |
| | ≤0.06 to ≥128 ⁷ | Enterobacterales | ≤1 | 2 | ≥4 |
| | ≤0.06 to ≥128 | <i>P. aeruginosa</i> | ≤2 | 4 | ≥8 |
| Meropenem-vaborbactam | ≤0.25 to ≥64 ⁸ | Enterobacterales | ≤4 | 8 | ≥16 |
| Piperacillin-tazobactam | ≤1 to ≥512 | <i>A. baumannii</i> complex | ≤16 | 32-64 | ≥128 |
| | ≤0.25 to ≥512 ⁹ | Enterobacterales | ≤8 | 16 | ≥32 |
| Tigecycline | ≤0.03 to ≥32 | Enterobacterales | ≤2 | 4 | ≥8 |
| Tobramycin | ≤0.06 to ≥64 ¹⁰ | Enterobacterales | ≤2 | 4 | ≥8 |
| | ≤0.125 to ≥64 | <i>P. aeruginosa</i> | ≤1 | 2 | ≥4 |
| Trimethoprim-sulfamethoxazole | ≤0.06 to ≥16 ¹¹ | Enterobacterales | ≤2 | - | ≥4 |

¹ AST Reportable range for Amikacin is ≤2 to ≥256 µg/mL for *C. koseri*, *E. coli*, *M. morgani* and *P. vulgaris*.

² AST Reportable range for Ampicillin-sulbactam is ≤2 to ≥128 µg/mL for *P. vulgaris*.

³ AST Reportable range for Aztreonam is ≤0.5 to ≥128 µg/mL for *C. freundii* complex, *E. coli* and *M. morgani*.

⁴ Interpretive category results for Cefepime/Enterobacterales are susceptible-dose dependent (SDD).

⁵ AST Reportable range for Ceftazidime is ≤0.5 to ≥128 µg/mL for *C. freundii* complex, *C. koseri*, *K. aerogenes* and *M. morgani*.

⁶ AST Reportable range for Ceftazidime-avibactam is ≤1 to ≥64 µg/mL for *E. coli*, *K. aerogenes*, *K. pneumoniae* group, *M. morgani* and *P. vulgaris*.

⁷ AST Reportable range for Meropenem is ≤0.125 to ≥64 µg/mL for *E. cloacae* complex, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae* group and *M. morgani*.

⁸ AST Reportable range for Meropenem-vaborbactam is ≤0.5 to ≥64 µg/mL for *M. morgani* and *P. vulgaris*.

⁹ AST Reportable range for Piperacillin-tazobactam is ≤1 to ≥512 µg/mL for *C. freundii* complex, *E. coli*, *K. aerogenes*, *K. pneumoniae* group and *M. morgani*.

¹⁰ AST Reportable range for Tobramycin is ≤0.125 to ≥64 µg/mL for *K. aerogenes*, *K. oxytoca*, *K. pneumoniae* group and *P. vulgaris*.

¹¹ AST Reportable range for Trimethoprim-sulfamethoxazole is ≤0.25 to ≥16 µg/mL for *C. freundii* complex, *C. koseri* and *M. morgani*.