

# A Pooled Performance Analysis of AS<sup>Tar</sup><sup>®</sup> Rapid Antimicrobial Susceptibility Testing System from *Pseudomonas aeruginosa* Positive and Contrived Blood Cultures

<sup>1</sup>Myers M, <sup>2</sup>Henry N, <sup>2</sup>Mutnal M, <sup>3</sup>Garrett E, <sup>3</sup>Bobenchik A, <sup>4</sup>Demessie G, <sup>4</sup>Yee R, <sup>5</sup>Cintrón M, <sup>6</sup>Marrero E, <sup>6</sup>Lauretta J, <sup>7</sup>Bao J, <sup>7</sup>Ponraj V, <sup>1</sup>Grawe J, <sup>1</sup>Gupta V\* \*Corresponding author: vikas.gupta@qlinea.com

<sup>1</sup>Q-linea, Uppsala, Sweden; <sup>2</sup>Baylor Scott and White Research Institute; <sup>3</sup>Penn State Hershey Medical Center, Hershey, PA; <sup>4</sup>George Washington University Hospital, Washington, DC; <sup>5</sup>Memorial Sloan Kettering Cancer Center, New York City, NY; <sup>6</sup>Lakeland Regional Health, Lakeland, FL; <sup>7</sup>Quest Diagnostics, Chantilly, VA

## Background

*Pseudomonas aeruginosa* (PsA) bloodstream infections are associated with high patient mortality and substantial healthcare burden. There is a serious need for rapid and accurate antimicrobial susceptibility testing (AST) when caring for patients with this type of bloodstream infection.

To assess performance of the rapid AST system AS<sup>Tar</sup><sup>®</sup> (Q-linea AB) (Figure 1) and the newly-FDA-cleared BC G-Kit v2 (Table 1), a pooled analysis was conducted using both clinical and contrived PsA-positive blood cultures (PBCs).

The analysis evaluated discrepancy testing between AS<sup>Tar</sup> and standard of care (SoC) methods, including: 1) repeat testing with SoC AST systems, and 2) adjudication with broth microdilution (BMD) as the reference method.



Figure 1. AS<sup>Tar</sup> Instrument and kit consumables (cartridge and AST disc) (not to scale)

## Methods

- Pooled AST results generated using AS<sup>Tar</sup> (7 US sites) with the AS<sup>Tar</sup> BC G-Kit panel (Investigational Use Only), were compared to SoC testing of PsA-positive blood cultures of overnight subcultures and contrived isolates (CDC AR Bank<sup>1</sup> and clinical sources). SoC testing was performed using Microscan WalkAway (4 sites) or Vitek 2 (3 sites).
- Performance was evaluated across 11 antipseudomonal antibiotics using 2025 FDA STIC interpretive criteria<sup>2</sup>: essential agreement (EA), categorical agreement (CA), very major errors (VME), and major errors (ME).
- Isolates with AST discrepancies underwent repeat testing using SoC methods and/or adjudication by BMD.
- Where available, BMD results from Sensititre (Thermo Fisher) were prioritized over BMD results from CDC AR Bank isolates.

## Conclusions

- Overall EA and CA agreement between AS<sup>Tar</sup> and SoC methods was > 90% for PsA-positive blood cultures, both pre- and post-adjudication.
- Post-adjudication VME and ME rates were within acceptable performance thresholds.
- Sixty percent (3 of 5) of evaluated VME discrepancies were resolved in favor of AS<sup>Tar</sup> after repeat SoC testing or adjudication by reference BMD.

## References

- CDC & FDA Antimicrobial Resistance Isolate Bank: <https://www.cdc.gov/arisolatebank/>
- U.S. Food and Drug Administration. Antibacterial Susceptibility Test Interpretive Criteria. Available at: <https://www.fda.gov/drugs/development-resources/antibacterial-susceptibility-test-interpretive-criteria>

## Results

A total of 71 PsA isolates were evaluated across 7 sites: 28.2% (20/71) were carbapenem-resistant, 21.1% (15/71) were multidrug-resistant, and 11.3% (8/71) were difficult-to-treat resistant.

High agreement was observed both pre- and post-adjudication:

- Overall EA of 96.3% (698/725) and 96.7% (701/725), pre- and post-adjudication
- Overall CA of 93.6% (644/688) and 94.3% (649/688), pre- and post-adjudication
- VME rates decreased from 3.7% (5/135) to 1.5% (2/133), following adjudication
- ME rates decreased from 1.0% (5/500) to 0.6% (3/500), following adjudication.

Overall, post-adjudication EA for all antibiotics was > 90%, and post-adjudication CA was > 90% for all antibiotics except aztreonam (86.9%) and cefepime (87.1%). Post-adjudication VME rates were < 3.0% for all antipseudomonal antibiotics, except for amikacin (33.3%) and tobramycin (12.5%) (Table 2).

Notably, 60% (3/5) of VMEs were resolved in favor of AS<sup>Tar</sup> following repeat SoC testing or BMD adjudication, including aztreonam, ceftazidime-avibactam, and piperacillin-tazobactam.

Table 1. April 2024 (v1), AS<sup>Tar</sup> BC G- FDA cleared panel.

Antimicrobial class	Antimicrobial agent	A. baumannii complex	C. freundii complex	C. koseri	E. cloacae complex	E. coli	K. aerogenes	K. oxytoca	K. pneumoniae group	M. morganii	P. mirabilis	P. vulgaris	S. marcescens	P. aeruginosa
Penicillin	Ampicillin													
β-lactam combination agents	Ampicillin-sulbactam <sup>1</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•
β-lactam combination agents	Ceftolozane-tazobactam <sup>2</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•
β-lactam combination agents	Ceftazidime-avibactam <sup>3</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•
β-lactam combination agents	Meropenem-vaborbactam <sup>4</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•
β-lactam combination agents	Piperacillin-tazobactam <sup>5</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Cefazolin		•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Cefepime	•	•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Cefotaxime	•	•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Ceftriaxone	•	•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Cefoxitin		•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Cefuroxime		•	•	•	•	•	•	•	•	•	•	•	•
Cephalosporin	Ceftazidime	•	•	•	•	•	•	•	•	•	•	•	•	•
Monobactam	Aztreonam		•	•	•	•	•	•	•	•	•	•	•	•
Carbapenem	Ertapenem		•	•	•	•	•	•	•	•	•	•	•	•
Carbapenem	Meropenem	•	•	•	•	•	•	•	•	•	•	•	•	•
Aminoglycoside	Gentamicin		•	•	•	•	•	•	•	•	•	•	•	•
Aminoglycoside	Tobramycin		•	•	•	•	•	•	•	•	•	•	•	•
Aminoglycoside	Amikacin	•	•	•	•	•	•	•	•	•	•	•	•	•
Tetracycline	Tigecycline		•	•	•	•	•	•	•	•	•	•	•	•
Fluoroquinolone	Ciprofloxacin		•	•	•	•	•	•	•	•	•	•	•	•
Fluoroquinolone	Levofloxacin		•	•	•	•	•	•	•	•	•	•	•	•
Miscellaneous	Trimethoprim-sulfamethoxazole <sup>5</sup>	•	•	•	•	•	•	•	•	•	•	•	•	•

<sup>1</sup>Ampicillin-sulbactam in the ratio 2:1

<sup>2</sup>For susceptibility testing purposes, the concentration of tazobactam is fixed at 4 μg/mL

<sup>3</sup>For susceptibility testing purposes, the concentration of avibactam is fixed at 4 μg/mL

<sup>4</sup>For susceptibility testing purposes, the concentration of vaborbactam is fixed at 8 μg/mL

<sup>5</sup>Trimethoprim-sulfamethoxazole in the ratio 1:19

Table 2. April 2026 (v2). Pre- and Post Adjudication Essential Agreement (EA), Categorical Agreement (CA), Very Major Errors (VME) and Major Errors (ME) Summary.

Antimicrobial	Pre-Adjudication				Post-Adjudication			
	EA % (n/N)	CA % (n/N)	VME % (n)	ME % (n)	EA % (n/N)	CA % (n/N)	VME % (n)	ME % (n)
Amikacin <sup>†</sup>	98.4% (62/63)	98.4% (62/63)	33.3% (1)	0% (0)	98.4% (62/63)	98.4% (62/63)	33.3% (1*)	0% (0)
Aztreonam	90.2% (55/61)	85.2% (52/61)	6.3% (1)	2.5% (1)	91.8% (56/61)	86.9% (53/61)	0% (0)	2.4% (1)
Cefepime	92.9% (65/70)	87.1% (61/70)	0% (0)	2.0% (1)	92.9% (65/70)	87.1% (61/70)	0% (0)	2.0% (1)
Ceftazidime <sup>†</sup>	92.2% (59/64)	96.9% (62/64)	0% (0)	0% (0)	92.2% (59/64)	96.9% (62/64)	0% (0)	0% (0)
Ceftazidime-avibactam	96.9% (63/65)	93.8% (61/65)	14.3% (1)	5.2% (3)	98.5% (64/65)	98.5% (64/65)	0% (0)	1.8% (1)
Ceftolozane-tazobactam	100% (64/64)	95.3% (61/64)	0% (0)	0% (0)	100% (64/64)	95.3% (61/64)	0% (0)	0% (0)
Ciprofloxacin	98.6% (69/70)	95.7% (67/70)	0% (0)	0% (0)	98.6% (69/70)	95.7% (67/70)	0% (0)	0% (0)
Levofloxacin	98.6% (69/70)	91.4% (64/70)	0% (0)	0% (0)	98.6% (69/70)	91.4% (64/70)	0% (0)	0% (0)
Meropenem	98.6% (69/70)	97.1% (68/70)	0% (0)	0% (0)	98.6% (69/70)	97.1% (68/70)	0% (0)	0% (0)
Piperacillin-tazobactam <sup>†</sup>	97.0% (65/67)	95.5% (64/67)	5.3% (1)	0% (0)	98.5% (66/67)	97.0% (65/67)	0% (0)	0% (0)
Tobramycin	95.1% (58/61)	91.7% (22/24)	12.5% (1)	0% (0)	95.1% (58/61)	91.7% (22/24)	12.5% (1*)	0% (0)
<b>Grand Total</b>	<b>96.3% (698/725)</b>	<b>93.6% (644/688)</b>	<b>3.7% (5)</b>	<b>1.0% (5)</b>	<b>96.7% (701/725)</b>	<b>94.3% (649/688)</b>	<b>1.5% (2*)</b>	<b>0.6% (3)</b>

\* Method of adjudication was repeat by comparator method and not adjudicated by BMD.

<sup>†</sup> Not FDA cleared